

Physicochemical characterization of manufactured nanomaterials

(TiO₂, SiO₂) used for genotoxicity testing


C. Guiot and O. Spalla, camille.guiot@cea.fr

DIRECTION DES SCIENCES DE LA MATIERE
Institut Rayonnement Matière de Saclay



European Joint Action 2010-2013

- Methodology to test genotoxicity of manufactured nanomaterials (MNs) (harmonized protocols)
- Relevant and reliable data for risk assessment by public health authorities



- Coordination (WP1)
- Dissemination (WP2)
- Evaluation (WP3)

❖ **WP4** (5 institutions – 4 countries)



❖ Detailed physico-chemical properties

❖ Characterization techniques

Size, shape, polydispersity, specific surface area, impurities, surface modification, suspension stability, etc.

XDR, Raman, TEM, SEM, **AFM**, BET, **SAXS**, ICP-MS/ AAS, MALDI-TOF, Dustiness, **DLS**, **zetametry**, redox, etc.

Small Angle X-ray Scattering (SAXS)

Radiation-matter interaction

Scattering vector
 $\vec{q} = \vec{k}_d - \vec{k}_i$
 $q = \frac{4\pi \sin\theta}{\lambda}$

SAXS measurements on TiO₂ powders

Specific surface area

Sample	Σ [m ²]	Specific surface area	Equivalent diameter for spheres
A2	2.78E+08	66 m ² /g	22 nm
R'4	2.22E+08	52 m ² /g	27 nm
RA5	1.99E+08	47 m ² /g	30 nm

SAXS measurements on SiO₂

Size
 ➤ Polydispersity
 ➤ Aggregation state
 ➤ Experimental data fitted by models of size distribution (nanostructure)

Stability and aggregation over time

Dynamic Light Scattering (DLS) - Zetametry

❖ Size measurement (hydrodynamic diameter) ❖ Zeta potential measurement (surface charge)

Stokes-Einstein
 $d_h = \frac{kT}{3\pi\eta D}$

Electrophoretic mobility
 $U_E = \frac{2\epsilon_0 \epsilon_r f(\kappa a)}{3\eta}$

Aggregation state

Protocol	Z-average (nm)	PDI
Protocol 1	380	0.303
Protocol 2	401	0.367
Protocol 3	685	0.336

Stability over pH range

Stability over time

Atomic Force Microscopy (AFM)

➤ **Z** ⇔ accurate measurement (nm), independent of the probe

➤ **x, y** ⇔ convolution of the probe diameter

Charged particles (positive) **Freshly cleaved mica sheet**

Spontaneous adsorption from a dilute suspension **Rinsing, ageing, drying**

Shape
 ➤ Aggregation state
 ➤ Height distribution

Statistics (size, polydispersity)



<http://www.nanogenotox.eu/>
<http://iramis.cea.fr/sis2m/lions>

* This communication arises from the Joint Action NANOGENOTOX which has received funding from the European Union, in the framework of the Health Programme. It reflects only the authors' views and the Executive Agency for Health and Consumers is not liable for any use that may be made of the information contained therein. The nanomaterials presented here come from the JRC repository and are also studied in the sponsorship program of OECD.

Dr. Camille GUIOT

Interdisciplinary Laboratory on Molecular Systems and Materials

DSM/IRAMIS/SIS2M/LIONS

Bat 125 point courrier 9

CEA Saclay

91191 Gif-sur-Yvette

France

Tel : +33 (0)1 69 08 43 59

Fax : +33 (0)1 69 08 66 40

camille.guiot@cea.fr